

UNICO® PowerSpin™ LX PRO Series Centrifuge

The C85608 PowerSpin™ LX PRO Centrifuge offers digital variable speed control from 300 to 4,000 RPM, providing flexibility for different sample types. With an 8-place rotor and a 99-minute digital timer, it can accommodate up to 8 x 15ml tubes, making it ideal for a wide range of laboratory applications. This centrifuge meets top safety and electromagnetic compatibility standards for reliable operation in laboratory environments.



LX Adapters

The LX PRO Centrifuge offers optional adapters for specialized sample types. The Microhematocrit adapter (set of 2) is designed for use with capillary tubes. The Short tube adapter (set of 2) is compatible with 3-5ml blood tubes (75mm in length), and the Pediatric or microtube adapter (pack of 2) accommodates 38-45mm length tubes, ideal for small-volume samples. These adapters provide additional flexibility, allowing the centrifuge to handle a wider variety of tubes and sample types.

Item No.	Description
C800-34	Microhematocrit adapter. Set of 2.
C800-06	Short tube adapter for 3-5ml (75mm) blood tubes. Set of 2.
C800-04	Microtube (38-45mm tube length), Set of 2.



8 Place Fixed Angle for 2 to 15ml tubes.



Short Tube Adapter (75mm Tube Length), pack of 6.



Microhematocrit adapter. Set of 2.



Short tube adapter for 3-5ml (75mm) blood tubes. Set of 2 pieces.



Adapter for pediatric or microtube (38-45mm tube length), pack of 2.

Specifications:

- Max. Rotation Speed: 4000rpm (300-4000rpm adjustable, 100rpm steps)
- Max. Relative Centrifugal Acceleration: 1860 xg
- Capacity: 15ml x 8 tubes
- Timing: 1 minute to 99 minutes, with continuous operation (Hold mode)
- Noise Level: ≤62dB(A) for quiet operation
- Drive Motor: Brushless DC motor for reliable, maintenance-free performance
- Safety Features: Includes door lock, overspeed protection, and status diagnosis system
- Power Supply: 100V-240V, 50Hz/60Hz
- Protection Class: Class I for safe operation
- EMC Compliance: Radiated interference and immunity certified to EN/IEC 61326-1 & FCC Class A

Item No.	Description
LX-C85608	CENTRIFUGE,LX, 8 PLACE FIXED ANGLE ROTOR, VARIABLE SPEED UP TO 4000RPM, 3 PRESET MODES.PLUG,110V220V50HZ60HZ

Shipping Dim: 22" x 16" x 18"

Net Weight: 20 lbs Shipping Weight: 22 lbs