


# High School Microscope

Model: M220FL



- MONOCULAR HEAD
- WF10X EYEPIECE
- 4X, 10X, 40X DIN ACHROMAT OBJECTIVES
- PLAIN STAGE WITH 0.65 CONDENSER
- FLUORESCENT ILLUMINATOR



The M220 Series of microscopes is designed for high school biology classrooms. These durable and reliable instruments are backed by a Limited Lifetime Warranty, ensuring long-lasting performance. Choose from a variety of models and features to suit your specific needs, including traditional tungsten illumination or energy-efficient cool-lite fluorescent options.

Available head options range from inclined monocular to dual "teaching" heads. All models are parfocal and parcentered, and all stages are drilled and tapped for optional mechanical stages. Select models are available with factory-installed mechanical stages. Each microscope includes an arm-mounted cord wrap, dust cover, and instruction manual.



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>ACCESSORIES</b>	<p><b>Miscellaneous Accessories</b></p> <p>M220-7007 AC ADAPTER FOR FREEDOM JR. SERIES</p> <p>M220-7003 BULB, FLUORESCENT, 5 WATT</p> <p>M220-5101 GRADUATED MECHANICAL STAGE, ATTACHABLE</p>
<b>SPECIFICATION</b>	<p><b>Head</b> Monocular</p> <p><b>Eyepiece</b> 10x WF</p> <p><b>Objectives</b> 4x, 10x, 40x</p> <p><b>Condenser</b> NA 0.65</p> <p><b>Diaphragm</b> Disc, 5-hole</p> <p><b>Stage</b> Plain, tapped</p> <p><b>Illumination</b> Fluorescent</p> <p><b>Bulb</b> 7 Watt</p> <p><b>Focusing</b> Coarse/Fine</p>
<b>SIZE AND WEIGHT</b>	<p><b>Actual Size</b> 8.5 in x 5 in x 11.5 in (216 mm x 127 mm x 292 mm)</p> <p><b>Shipping Size</b> 12 in x 8 in x 16 in (305 mm x 203 mm x 406 mm)</p> <p><b>Weight</b> 12 lbs (3.6 kg)</p> <p><b>Shipping Weight</b> 13 lbs (4.1 kg)</p>
<b>WARRANTY</b>	1 years

