


PowerSpin™ PX Centrifuge

Model: C827



- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- POWER INPUT 15VDC



The C827 PXV 12V Variable Speed Centrifuge is a versatile and high-performance solution for mobile and field applications. Designed for 12V DC operation, it includes a vehicle power adapter (car adapter) for convenient use in remote locations, on-site testing, or mobile laboratories. This centrifuge offers variable speed control, allowing you to adjust the speed according to your specific sample processing needs. It features a 6-place rotor, accommodating up to six tubes for efficient sample separation. The dual lid lock ensures enhanced safety by securely locking the lid in place during operation, while the built-in leveling system provides stability even on uneven surfaces, ensuring consistent performance. With the C827 PXV, you can count on reliable and efficient centrifugation wherever you go, combining convenience and performance for a wide range of applications.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories (Sold Separately)</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p> <p>C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement</p> <p>C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord</p>
SPECIFICATION	<p>Number of Tubes 6</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Variable</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3800</p> <p>Maximum RCF(g) 1320</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p>



	Weight	18 lbs
	Shipping Weight	19 lbs
WARRANTY		1 years

