

PowerSpin™ PX Centrifuge

Model: C828H



- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- FIXED SPEED 3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 15VDC



The C828 PowerSpin™ PX Centrifuge is an advanced, high-performance portable centrifuge designed to meet the demands of diverse laboratory applications. Equipped with an eight-place rotor, it is compatible with standard 2-15 ml test tubes, as well as 2 ml and 4-7.5 ml tubes when used with optional adapters. This versatility makes it ideal for handling a variety of sample sizes in clinical, research, and industrial settings. Designed for 12V DC operation, it includes a vehicle power adapter (car adapter) for convenient use in remote locations, on-site testing, or mobile laboratories. This centrifuge offers variable speed control, allowing you to adjust the speed according to your specific sample processing needs. It features a 8-place rotor, accommodating up to eight tubes for efficient sample separation. The dual lid lock ensures enhanced safety by securely locking the lid in place during operation, while the built-in leveling system provides stability even on uneven surfaces, ensuring consistent performance. With the C827 PXV, you can count on reliable and efficient centrifugation wherever you go, combining convenience and performance for a wide range of applications.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories (Sold Separately)</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p> <p>C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement</p> <p>C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord</p>
SPECIFICATION	<p>Number of Tubes 8</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Fixed</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3800</p> <p>Maximum RCF(g) 1320</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p>



	Weight	18 lbs
	Shipping Weight	19 lbs
WARRANTY		1 years

