

PowerSpin™ PX Centrifuge

Model: C829H



- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 15VDC



The C829 PX 12V Centrifuge is a robust and portable centrifuge designed to deliver reliable performance in mobile and field settings. Powered by a 12V DC vehicle power adapter (included), this centrifuge allows for easy operation directly from your vehicle's power source, making it perfect for remote locations, fieldwork, or mobile laboratories. With an 8-place rotor, the C829 PX accommodates up to eight tubes, enabling efficient processing of multiple samples at once. For added safety, it features a dual lid lock to ensure that the lid remains securely closed during operation. The built-in leveling system provides stability, ensuring consistent and reliable performance even on uneven surfaces. Whether you are conducting research in the field or need a centrifuge for on-the-go testing, the C829 PX offers the ideal solution, combining convenience, safety, and high performance in one portable unit.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p> <p>C826-06 Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement</p> <p>C826-12 Porta-Spin Ac Adapter, 110 Volt, USA Power Cord</p>
SPECIFICATION	<p>Number of Tubes 8</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Variable</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3800</p> <p>Maximum RCF(g) 1647</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p>



	Weight	18 lbs
	Shipping Weight	19 lbs
WARRANTY		1 years

