

PowerSpin™ HX Centrifuge

Model: C8304



- POWERSPIN™ HX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 12 PLACE FIXED ANGLE ROTOR
- 30 MIN. DIGITAL TIMER
- 6 X10ML CAPACITY
- POWER INPUT 110V



The UNICO PowerSpin™ HX Series benchtop centrifuges provide an economical, reliable, and exceptionally quiet centrifugation solution, available with either horizontal or 12-place fixed angle rotors. Utilizing a horizontal rotor ensures a clean and complete separation of cellular components, yielding superior serum or plasma specimens. Specifically, when employing a horizontal rotor with Serum Separator Tubes (SSTs), a tight, uniform horizontal seal is created, effectively minimizing the risk of red cell leakage and potential specimen deterioration, a critical factor for specimens sent for external testing. The PowerSpin™ HX model incorporates an electronically controlled safety lock, preventing lid opening during rotor motion, along with a digital timer and LCD display for precise spin cycle control. Its powerful, brushless motor ensures years of dependable, maintenance-free operation, and the stainless steel guard bowl facilitates easy cleaning and maintenance. The PowerSpin™ HXV model features a Zero RPM safety lock, which also prevents the lid from being opened while the rotor is in motion, and it also contains a digital timer and LCD display for precise spin cycle control.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Accessories (Sold Separately) C8300-01 6 place horizontal rotor for HXV C8300 centrifuge C8300-02 8 place horizontal rotor for HXV C8300 centrifuge C8300-03 12 place angled rotor for HXV C8300 centrifuge
SPECIFICATION	Speed (RPM) 3500 Maximum RCF (g) 1777g Timer Digital Nominal Capacity 12 x 15 ml Rotor Fixed Angle Lock ZERO RPM Digital Motor AC Maximum Test Tube Size 16 x 110 mm
WARRANTY	1 years

