


PowerSpin™ HX Centrifuge

Model: C8306



- POWERSPIN™ HXV CENTRIFUGE
- VARIABLE SPEED 500-3500 RPM
- 6 PLACE HORIZONTAL ROTOR
- 99 MIN. DIGITAL TIMER
- 6 X10ML CAPACITY
- POWER INPUT 110V



The UNICO® PowerSpin™ HX Series Centrifuge (Model # C8306) is a highly reliable and versatile laboratory tool, designed to meet the demands of clinical, industrial, and biological labs. It features a 6-place swing-out rotor, perfect for 3-15 ml tubes, ensuring efficient separation of blood and other materials. With easy-to-use push button keypad controls, the centrifuge can be operated with gloves on, making it convenient for busy lab environments.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Accessories (Sold Separately) C8300-01 6 place horizontal rotor for HXV C8300 centrifuge C8300-02 8 place horizontal rotor for HXV C8300 centrifuge C8300-03 12 place angled rotor for HXV C8300 centrifuge
SPECIFICATION	Speed (RPM) 500-3500 Maximum RCF (g) 1777g Timer Digital Nominal Capacity 6 x 15 ml Rotor Horizontal Lock ZERO RPM Digital Motor AC Maximum Test Tube Size 16 x 110 mm
WARRANTY	1 years

