

PowerSpin™ DX Centrifuge

Model: C8706



- POWERSPIN™ DX CENTRIFUGE
- VARIABLE SPEED 500-4000 RPM
- PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL
- 6 PLACE HORIZONTAL ROTOR
- 1:00-99:99 MIN DIGITAL TIMER
- CAPACITY 6X15 ML TUBES
- POWER INPUT 110V



The UNICO PowerSpin™ DX Centrifuge is engineered to meet the demanding needs of busy veterinary clinics. Featuring digital variable speed control (500-4000 RPM), an electronic timer, and a user-friendly interface with preset programs for common applications, the C8706 ensures efficient and reliable sample processing. Its compact design and quiet operation minimize disruptions in your clinic, while the maintenance-free brushless DC motor ensures long-lasting performance.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories (Sold Separately)</p> <p>C8600-01 6 Place Horizontal Rotor, 2-15ml for DX Centrifuge</p> <p>C8600-02 4 Place Horizontal rotor, 50ml DX Centrifuge</p> <p>C8600-03M 24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)</p>
SPECIFICATION	<p>Speed (RPM) 500-4000rpm</p> <p>Maximum RCF (g) 1848g</p> <p>Rotor 6 x 15ml Horizontal</p> <p>Tube Size 2-15ml</p> <p>Motor DC</p> <p>Power Input 110V</p>
SIZE AND WEIGHT	<p>Size 13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width</p> <p>Shipping Size 19"X18"X15"</p> <p>Weight 44 LBS</p> <p>Shipping Weight 44 LBS</p>
WARRANTY	1 years

