


PowerSpin™ BX Centrifuge

Model: C881



- POWERSPIN™ BX CENTRIFUGE
- VARIABLE SPEED 1000-5000 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10ML OR 8X12ML TUBES
- POWER INPUT 110V



Boost your lab's efficiency with the C881 BX High-Speed Centrifuge. This robust benchtop centrifuge boasts an 8-place rotor, accommodating 8x10ml tubes, and achieves speeds up to 5000 RPM (up to 11,000 RPM with compatible rotors). The oversized LCD display provides clear, real-time data, while the quiet, maintenance-free brushless DC motor ensures reliable operation. A unique feature is the integrated reader on the microhematocrit rotor cover, eliminating the need for separate reading devices. Operates on 110V.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories (Sold Separately)</p> <p>C880-08 Fixed Angle Rotor, 8 place</p> <p>C880-02 Microhematocrit Rotor, 24 place</p> <p>C880-03 Microtube rotor, 24 place</p> <p>C880-10 BX Rotor Rack</p> <p>C881-03 Short tube adapter, for 13x75mm tubes</p> <p>C881-04 Micro / Pediatric tube adapter for 0.5-2.0ml tubes</p>
SPECIFICATION	<p>Rotor Fixed Angle, 8 place rotor for 10ml SST Tubes</p> <p>Speed Range 1,000 to 5,000 RPM</p> <p>Maximum RCF(g) 2571g</p> <p>Timer Range 0-99 minutes LCD</p> <p>Max. Rotor Capacity 8 x 10ml or 8x12ml</p> <p>Dimension 270mm x 320mm x 220mm</p>
SIZE AND WEIGHT	<p>Size 10.5 in x 12.5 in x 9 in (270 mm x 320 mm x 220 mm)</p> <p>Shipping Size 19 in x 15 in x 13 in</p> <p>Weight 15.5 lbs (C881 & C882); 16.5 lbs (C883)</p> <p>Shipping Weight 19 lbs (C881 and C882); 20 lbs (C883)</p>
WARRANTY	1 years

