

PowerSpin™ BX Centrifuge

Model: C884



- POWERSPIN™ BX CENTRIFUGE
- VARIABLE SPEED 1000-11000 RPM
- 8 PLACE FIXED ANGLE AND 24PLACE MICROHEMATOCRIT ROTOR
- 99 MIN DIGITAL TIMER
- POWER INPUT 110V



The BX Centrifuge is a versatile benchtop centrifuge designed for a range of laboratory applications. This particular model offers exceptional flexibility with two included rotors: a larger 8-place rotor designed for 10mL tubes (8 x 10mL capacity) ideal for larger volume sample processing, and a 24-place microhematocrit rotor for efficient and accurate hematocrit determinations. This dual-rotor capability eliminates the need for purchasing separate centrifuges, saving valuable lab space and budget. The BX Centrifuge is powered by a standard 110V power supply, making it readily adaptable to most laboratory environments. The built in oversized LCD Display provides real time data to the user. With the push of a button you can speed up or slow down the rotor within 1,000 - 11,000 RPM range (Max 5000 RPM with the Blood Tube Rotor). Maintenance free brushless DC motor is powerful and super quiet. The cover of the Microhematocrit rotor incorporates a reader, allowing for instant reading and eliminating the need for additional specimen handling.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories (Sold Separately)</p> <p>C880-08 Fixed Angle Rotor, 8 place</p> <p>C880-02 Microhematocrit Rotor, 24 place</p> <p>C880-03 Microtube rotor, 24 place</p> <p>C880-10 BX Rotor Rack</p> <p>C881-03 Short tube adapter, for 13x75mm tubes</p> <p>C881-04 Micro / Pediatric tube adapter for 0.5-2.0ml tubes</p>
SPECIFICATION	<p>Rotor Microcentrifuge, 24 place</p> <p>Speed Range 1000 - 11000</p> <p>Maximum RCF(g) 11093</p> <p>Nominal Capacity 24 x 0.5-2.0 ml</p> <p>Maximum Test Tube Size 5 x 10mm</p>
SIZE AND WEIGHT	<p>Size 10.5 in x 12.5 in x 9 in (270 mm x 320 mm x 220 mm)</p> <p>Shipping Size 14 in x 15 in x 12 in (350 mm x 375 mm x 300 mm)</p> <p>Weight 15.5 lbs (C881 & C882); 16.5 lbs (C883)</p> <p>Shipping Weight 19 lbs (C881 and C882) ; 20 lbs (C883)</p>
WARRANTY	1 years

