

PowerSpin™ BX Centrifuge

Model: C885



- POWERSPIN™ BX CENTRIFUGE
- VARIABLE SPEED 1000-13000 RPM WITH 8 PLACE FIXED ANGLE, 24 PLACE MICROHEMATOCRIT AND 24 PLACE MICROCENTRIFUGE ROTORS
- 99 MIN DIGITAL TIMER
- POWER INPUT 110 V



The BX Centrifuge is a highly versatile and efficient benchtop centrifuge designed to meet the diverse needs of modern laboratories. This model comes equipped with three interchangeable rotors for maximum flexibility: a 24-place microcentrifuge rotor (24 x 1.5-2mL), a 24-place microhematocrit rotor, and an 8-place rotor for 10mL tubes (8 x 10mL). A key feature is the automatic rotor detection system, which simplifies operation and ensures proper settings for each rotor. The included rotor changing tool makes switching between rotors quick and easy. This multi-rotor capability eliminates the need for multiple centrifuges, saving valuable bench space and budget. The BX Centrifuge operates on standard 110V power. The built in oversized LCD Display provides real time data to the user. With the push of a button you can speed up or slow down the rotor within 1,000 - 11,000 RPM range (Max 5000 RPM with the Blood Tube Rotor). Maintenance free brushless DC motor is powerful and super quiet. The cover of the Microhematocrit rotor incorporates a reader, allowing for instant reading and eliminating the need for additional specimen handling.



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories (Sold Separately) C880-08 Fixed Angle Rotor, 8 place C880-02 Microhematocrit Rotor, 24 place C880-03 Microtube rotor, 24 place C880-10 BX Rotor Rack C881-03 Short tube adapter, for 13x75mm tubes C881-04 Micro / Pediatric tube adapter for 0.5-2.0ml tubes	
SPECIFICATION	Power Input	110V
SIZE AND WEIGHT	Size	10.5 in x 12.5 in x 9 in (270 mm x 320 mm x 220 mm)
	Shipping Size	14 in x 15 in x 12 in (350 mm x 375 mm x 300 mm)
	Weight	15.5 lbs (C881 & C882); 16.5 lbs (C883)
	Shipping Weight	19 lbs (C881 and C882) ; 20 lbs (C883)
WARRANTY	1 years	

