


# PowerSpin™ MH Centrifuge


Model: CMH30



- POWERSPIN™ MH MICRO-HEMATOCRIT CENTRIFUGE
- FIXED SPEED 12000 RPM
- 24 PLACE MICROHEMATOCRIT ROTOR
- 5 MIN DIGITAL TIMER
- POWER INPUT 110V
- MICROHEMATOCRIT READER INCLUDED.



The CMH30 Heavy Duty 24-Place Micro Hematocrit Centrifuge is a reliable and efficient solution for determining hematocrit levels in blood samples. This 110V model is built for demanding laboratory environments and features a robust design capable of handling high volumes. It accommodates up to 24 micro hematocrit capillary tubes, allowing for simultaneous processing and increased throughput. The integrated 5-minute timer ensures consistent and accurate spin cycles, while the convenient braking mechanism allows for quick deceleration and faster turnaround times. This package includes a built-in reader for easy and precise hematocrit determination, eliminating the need for separate equipment. For added convenience and extended use, an extra gasket is also included.



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>ACCESSORIES</b>	<b>Accessories</b> CMH30-Reader Microhematocrit reader CMH30-03 Gasket Replacement For CMH30 Centrifuge
<b>SPECIFICATION</b>	<b>Speed</b> 12000 RPM <b>RCF</b> 14,837g <b>Rotor</b> 24 place <b>Power Input</b> 110V
<b>SIZE AND WEIGHT</b>	<b>Size</b> 9in x 10 in x 9 in <b>Shipping Size</b> 13 in x 13 in x 13 in <b>Weight</b> 18 lbs <b>Shipping Weight</b> 19.4 lbs
<b>WARRANTY</b>	1 years

