


# Differential Counter

Model: L-BC3



- UNICO DIFFERENTIAL COUNTER
- COMPACT DESIGN
- FEATURING 2-KEY AND RESET KNOB



This manual differential cell counter is a reliable and efficient tool for performing manual blood cell differentials. Featuring a two-key design, it simplifies the counting process and helps prevent errors. The integrated totalizer window displays the total count, while the audible signal bell confirms each count, allowing the user to focus on the slide. Dual reset knobs provide convenient clearing of individual and total counts. This counter is ideal for use in clinical laboratories and research settings.



**ADDITIONAL INFORMATION**

<b>TOP SELECTION</b>	No
<b>WARRANTY</b>	1 years

