

Differential Counter

Model: L-BC6



- UNICO DIFFERENTIAL COUNTER
- DURABLE ABS PLASTIC HOUSING
- 5-KEY WITH TOTALIZER WINDOW
- SIGNAL BELL AT 100 COUNT
- DUAL RESET KNOBS



The L-BC6 Differential Counter is a valuable tool for performing manual blood cell counts and other differential counting tasks in laboratory settings. This robust counter features five individual keys, each corresponding to a specific cell type (e.g., neutrophils, lymphocytes, monocytes, eosinophils, basophils), allowing for quick and accurate tallying. A totalizer window displays the overall count, while an audible signal bell rings every 100 cells counted, simplifying the counting process and reducing the risk of errors. Dual reset knobs allow for independent resetting of both the individual cell type counts and the totalizer, providing flexibility and convenience during use. The L-BC6 is a reliable and user-friendly device for accurate manual differential cell counting.



ADDITIONAL INFORMATION

TOP SELECTION	No
WARRANTY	1 years

