

# Differential Counter

Model: L-BC9



- UNICO DIFFERENTIAL COUNTER
- METAL HOUSING
- 8-KEY WITH TOTALIZER WINDOW
- WBC MATURATION PICTO-STRIP
- SIGNAL BELL AT 100 COUNT
- DUAL RESET KNOBS



The L-BC9 Differential Counter is a valuable tool for performing manual blood cell counts and other differential counting tasks, offering enhanced functionality with eight individual keys and a helpful maturation chart. This robust counter features eight individual keys, each corresponding to a specific cell type or stage of maturation in white blood cells, allowing for more detailed and accurate tallying. A totalizer window displays the overall count, while an audible signal bell rings every 100 cells counted, simplifying the counting process and reducing the risk of errors. A full white blood cell maturation strip chart is integrated into the device, providing a quick visual reference for cell identification. Dual reset knobs allow for independent resetting of both the individual cell type counts and the totalizer, providing flexibility and convenience during use. The L-BC9 is a reliable and user-friendly device for accurate and detailed manual differential cell counting, particularly useful in hematology and clinical laboratory settings.



**ADDITIONAL INFORMATION**

<b>TOP SELECTION</b>	No
<b>WARRANTY</b>	1 years

