

Tube Rocker

Model: L-TTR200D



- UNICO TEST TUBE ROCKER
- UP TO 30 TUBES/VIALS CAPACITY
- EASY OPERATION TOP PAUSE BUTTON
- FIXED SPEED OF 25 ROCKING CYCLES PER MINUTE
- POWER INPUT 110 V



This Double Decker Blood Tube Rocker significantly increases sample throughput with its 30-place capacity. Designed for efficient and gentle mixing, it accommodates a large number of blood tubes or vials simultaneously. The unit operates at a consistent 25 rocking cycles per minute (RPM) on standard 110V power, ensuring thorough mixing without agitation or damage to samples. Its double-decker design maximizes capacity while minimizing footprint, making it ideal for busy laboratories.



ADDITIONAL INFORMATION

TOP SELECTION	No
SPECIFICATION	Capacity 30 tubes Speed 25cpm Voltage 110V
SIZE AND WEIGHT	Size 8.6 in x 4.7 in x 5.0 in Shipping Size 10.6 in x 6.7 in x 6.9 in Weight 4.0 lbs (1.8 kg) Shipping Weight 5.0 lbs (2.3 kg)
WARRANTY	1 years

