

Vortex Mixer

Model: L-VM1000



- UNICO VORTEX MIXER
- SINGLE TUBE HEAD
- 3-POSITION SWITCH FOR CONSTANT OR "TOUCH-ON" OPERATION
- VARIABLE SPEED UP TO 3000RPM
- POWER INPUT 120V



The UNICO single tube vortex mixer is rugged and durable, designed for your most demanding lab applications. It features a dual-mode switch that lets you vortex continuously or use touch mode — simply press your tube against the soft rubber mixing cup to start intermittent mixing. The built-in speed control allows precise adjustment of the vortexing speed to suit your requirements.



ADDITIONAL INFORMATION

TOP SELECTION	No
SPECIFICATION	Tube Head Single Switch 3-Position for Constant or "Touch-On" Operation Variable Speed Up to 3000 RPM Power Input 120V
SIZE AND WEIGHT	Size 10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm) Shipping Size 12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm) Weight 9.8 lbs (4.4 kg) Shipping Weight 10.8 lbs (4.9 kg)
WARRANTY	1 years

