

# Vortex Mixer

Model: L-VM2000



- UNICO VORTEX MIXER
- 3.75 INCHES VIBRATION PAD
- 3-POSITION SWITCH FOR CONSTANT OR "TOUCH-ON" OPERATION
- VARIABLE SPEED UP TO 3000RPM
- POWER INPUT 120V



UNICO economical vortex mixers are rugged and sturdy for your most demanding applications. Dual mode switch lets you vortex continuously or intermittently using the touch mode - simply press your tube to the soft rubber mixing cup. Both model feature speed control to adjust the vortexing velocity to your requirements.



**ADDITIONAL INFORMATION**

<b>TOP SELECTION</b>	No
<b>SIZE AND WEIGHT</b>	<b>Size</b> 10.4 in x 9.2 in x 6.0 in (263 x 234 x 153 mm) <b>Shipping Size</b> 12.2 in x 11.2 in x 7.4 in (311 x 285 x 187 mm) <b>Weight</b> 9.8 lbs (4.4 kg) <b>Shipping Weight</b> 10.8 lbs (4.9 kg)
<b>WARRANTY</b>	1 years

