


Clinical Incubator

Model: L-CU100



- UNICO CLINICAL INCUBATOR
- 10L CHAMBER CAPACITY
- REMOVABLE SHELVE
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 110 VOLT



The L-CU100 Laboratory Incubator is a compact and reliable unit designed for a variety of incubation applications in research, clinical, and educational settings. With a 10-liter capacity, it's ideal for smaller-scale experiments and applications where space is limited. The incubator features precise digital temperature control and a clear LED display, allowing for accurate temperature setting and monitoring. Its 110V power input makes it suitable for standard North American electrical systems. The L-CU100 provides a stable and controlled environment for cell cultures, microbiological studies, and other temperature-sensitive applications.



ADDITIONAL INFORMATION

| | |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| TOP SELECTION | No |
| SPECIFICATION | <p>Chamber Capacity 10 Liters (0.35 cubic feet)</p> <p>Chamber Size Internal Dimensions 11.2"W x 7.8"D x 8.6"H</p> <p>Calibration Digital control & thermometer access port</p> <p>Temperature Range Ambient to 70°C (Ambient to 155° F)</p> <p>Shelves Removable stainless steel shelf</p> <p>Temperature Control Digital Control</p> <p>Display Digital Display</p> <p>Power/Consumption 1.4 Amps / 150 Watts</p> <p>Heating Convection</p> <p>Stability Will hold constant temperature +/- 0.5° C</p> |
| WARRANTY | 1 years |

