


# Clinical Laboratory Incubator

Model: L-CU200



- UNICO CLINICAL LABORATORY INCUBATOR
- 20L CHAMBER CAPACITY
- 2 REMOVABLE SHALVES
- TEMPERATURE AMBIENT TO 70° C
- POWER INPUT 110 VOLT



The L-CU200 Laboratory Incubator is a reliable and versatile unit designed for a variety of incubation applications in research, clinical, and educational settings. With a 20-liter capacity, it offers increased space for larger samples or higher throughput compared to smaller models. The incubator features precise digital temperature control and a clear LED display, allowing for accurate temperature setting and monitoring. Its 110V power input makes it suitable for standard North American electrical systems. The L-CU200 provides a stable and controlled environment for cell cultures, microbiological studies, and other temperature-sensitive applications requiring a larger incubation volume.



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>SPECIFICATION</b>	<p><b>Chamber Capacity</b> 20 Liters (0.7 cubic feet)</p> <p><b>Chamber Size</b> Internal Dimensions 12.3"W x 12.1"D x 9"H</p> <p><b>Calibration</b> Digital control &amp; thermometer access port</p> <p><b>Temperature Range</b> Ambient to 70°C (Ambient to 155° F)</p> <p><b>Shelves</b> Two removable shelves, perforated for air flow</p> <p><b>Temperature Control</b> Digital Control</p> <p><b>Display</b> Digital Display</p> <p><b>Power/Consumption</b> 1.4 Amps / 150 Watts</p> <p><b>Heating</b> Convection</p> <p><b>Stability</b> Will hold constant temperature +/- 0.5° C</p> <p><b>Power input</b> Available in 110V and 220V</p>
<b>WARRANTY</b>	1 years

