

# Med/Vet Practice Microscope

Model: G380SP-LED



- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X, 10X, 40XR, 100XR DIN SEMI-PLAN ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- LED ILLUMINATION



The G380SP-LED Binocular Microscope combines precision optics with a sleek, ergonomic design, making it an ideal choice for professionals and researchers in various laboratory settings. Equipped with a Seidentopf Binocular Head and 10x widefield eyepieces, this microscope ensures clear, comfortable viewing, even for users who wear glasses.

This model features a high-quality Semi-Plan Objective Set, which includes 4x, 10x, 40x, and 100x objectives. The semi-plan objectives provide enhanced image flatness and resolution across a wide range of magnifications, offering excellent performance for both routine and detailed observations.

The G380SP-LED comes with a mechanical stage, allowing for precise specimen manipulation with smooth, low-position coaxial controls. This feature ensures accurate movement and positioning of slides, providing the user with a stable and consistent experience.

For optimal lighting, the microscope includes built-in LED illumination with an intensity adjustment feature. This allows you to easily control the light levels for different specimens and applications, ensuring clear, bright images in any environment.

Whether you're conducting routine inspections or in-depth research, the G380SP-LED Binocular Microscope offers high-quality imaging, precise control, and flexible features to support your work.



## ADDITIONAL INFORMATION

TOP SELECTION	No	
<b>ACCESSORIES</b>	<b>Bodies</b>	
	G380-3001	MONOCULAR HEAD, 45° INCLINED FOR G380
	G380-3002	DUAL VIEWING HEAD WITHOUT EYEPIECE FOR G380 SERIES
	G380-3004	SEIDENTOPF BINOCULAR HEAD FOR G380 SERIES
	G380-3005	SEIDENTOPF TRINOCULAR HEAD FOR G380 SERIES
	<b>Optics</b>	
	G380-1402	HIGH EYEPOINT 10X WIDEFIELD EYEPIECE, F.N. 18, EACH
	G380-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10 FOR G380 SERIES
	G380-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25 FOR G380 SERIES
	G380-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS FOR G380 SERIES
	G380-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS FOR G380 SERIES
	G380-2301	4X DIN PLAN OBJECTIVE, N.A. 0.10
	G380-2302	10X DIN PLAN OBJECTIVE, N.A. 0.25
	G380-2305	25X DIN PLAN OBJECTIVE N.A. 0.40
	G380-2303	40X DIN PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	G380-2304	100X DIN PLAN OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS
	G380-2601	4X INFINITY ACHROMAT OBJECTIVE, NA 0.10
	G380-2602	10X INFINITY ACHROMAT OBJECTIVE, NA 0.25
	G380-2603	40X INFINITY ACHROMAT OBJECTIVE, NA 0.65
	G380-2604	100X OIL IMMERSION INFINITY ACHROMAT OBJECTIVE, NA 1.25
	G380-2631	4X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.10
	G380-2632	10X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.25



	G380-2633	40X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.65
	G380-2634	100X OIL IMMERSION INFINITY PLAN ACHROMAT OBJECTIVE, NA 1.25
	<b>Miscellaneous Accessories</b>	
	G380-5111	GRADUATED MECHANICAL STAGE (COAXIAL LOW POSITION CONTROLS) FOR G380 SERIES
	G380-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT FILTER HOLDER FOR G300 SERIES
	G380-7001	LED ILLUMINATOR PACK, SNAP-IN, INCLUDED 3 AA BATTERIES, FOR G380 SERIES ONLY
	G380-7002	MIRROR ASSEMBLY, SNAP-IN FOR G380 SERIES ONLY
	G380-7010	FIELD DIAGRAM KOEHLER SET FOR G380, G390 SERIES MICROSCOPES
	G380-8010	PHOTO LENS AND PHOTO ADAPTER WITH UNIVERSAL THREAD FOR G380 SERIES (REQ. CAMERA AND T-MOUNT)
	G380-8011	VERTICAL VIEWING ADAPTER (TO BE INSTALLED ON TRINOCULAR HEAD. REQUIRES EYEPIECE)
	G380-8031	PHASE KIT: INCLUDES 40XR SEMI-PLAN PHASE OBJECTIVE AND SLIDE-IN ANNULUS RING FOR G380 SERIES
	G380-8034	PHASE KIT: INCLUDES 10X SEMI-PLAN PHASE OBJECTIVE AND SLIDE-IN ANNULUS RING FOR G380 SERIES
	G380-8032	SLIDE-IN 40X PHASE ANNULUS RING FOR G380 SERIES
	G380-8907	HARD CARRYING CASE, PADDED WITH LOCK & KEY (2)
	G380-8001	DUST COVER FOR G380 SERIES
	G380-8906	DARKFIELD STOP FOR G380 SERIES
	G380-9008	BACKPACK™ FOR G380, G390 AND G490 SERIES MICROSCOPES, ACRYLIC ACCESORY ATTACHMENT
<b>SPECIFICATION</b>	<b>Head</b>	Binocular
	<b>Eyepiece</b>	10X WF Hi-EP
	<b>Objectives</b>	Semi-Plan 4X/10X/40X/100X
	<b>Condenser</b>	NA1.25
	<b>Diaphragm</b>	Iris
	<b>Stage</b>	Mechanical



	<b>Illumination</b>	LED
	<b>Bulb</b>	3 Watt LED
	<b>Focusing</b>	Coaxial
<b>SIZE AND WEIGHT</b>	<b>Shipping Size</b>	17 in x 10.5 in x 17.5 in
	<b>Weight</b>	17 lbs (7.7 kg)
	<b>Shipping Weight</b>	18 lbs (8.2 kg)
<b>WARRANTY</b>		1 years

