

Med/Vet Practice Microscope

Model: G383PL-LED



- SEIDENTOPF TRINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X, 10X, 40XR, 100XR DIN PLAN ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- LED ILLUMINATION



The G383PL-LED Trinocular Microscope combines advanced optical technology and user-friendly features, making it ideal for educational, laboratory, and research applications. Equipped with a Seidentopf trinocular head, this microscope offers simultaneous viewing for two users, with the additional option to capture high-quality images or video through the trinocular port. Perfect for collaborative learning and documentation, the G383PL-LED enhances both individual and group microscopy experiences.

This model includes 10x widefield eyepieces, ensuring clear, comfortable viewing with sharp, crisp images. The Plan Achromat objective set (4x, 10x, 40x, and 100x) delivers high-resolution, flat-field images, making it ideal for high-precision work where clarity and accuracy are essential. The Plan Achromat objectives minimize chromatic aberration and deliver excellent contrast across all magnifications, providing superior image quality for a variety of specimens.

The G383PL-LED also features a mechanical stage for precise specimen handling and manipulation. The Abbe condenser with an iris diaphragm allows for optimal light control, ensuring proper contrast and illumination for different types of specimens. This setup is perfect for a wide range of observation types, from standard brightfield to more detailed imaging.

For exceptional illumination, the G383PL-LED includes built-in LED lighting with intensity adjustment, ensuring uniform lighting across the specimen while minimizing heat generation. This energy-efficient LED system offers bright and consistent illumination, perfect for prolonged usage with minimal eye strain.

With its advanced optical features, ergonomic design, and versatile illumination options, the G383PL-LED Trinocular Microscope is a high-performance instrument for professionals and students alike, providing exceptional imaging for routine inspections, research, and educational use.





LIFETIME
LIMITED WARRANTY



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Bodies	
	G380-3001	MONOCULAR HEAD, 45° INCLINED FOR G380
	G380-3002	DUAL VIEWING HEAD WITHOUT EYEPIECE FOR G380 SERIES
	G380-3004	SEIDENTOPF BINOCULAR HEAD FOR G380 SERIES
	G380-3005	SEIDENTOPF TRINOCULAR HEAD FOR G380 SERIES
	Optics	
	G380-1402	HIGH EYEPOINT 10X WIDEFIELD EYEPIECE, F.N. 18, EACH
	G380-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10 FOR G380 SERIES
	G380-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25 FOR G380 SERIES
	G380-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS FOR G380 SERIES
	G380-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS FOR G380 SERIES
	G380-2301	4X DIN PLAN OBJECTIVE, N.A. 0.10
	G380-2302	10X DIN PLAN OBJECTIVE, N.A. 0.25
	G380-2305	25X DIN PLAN OBJECTIVE N.A. 0.40
	G380-2303	40X DIN PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	G380-2304	100X DIN PLAN OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS
	G380-2601	4X INFINITY ACHROMAT OBJECTIVE, NA 0.10
	G380-2602	10X INFINITY ACHROMAT OBJECTIVE, NA 0.25
	G380-2603	40X INFINITY ACHROMAT OBJECTIVE, NA 0.65
	G380-2604	100X OIL IMMERSION INFINITY ACHROMAT OBJECTIVE, NA 1.25
	G380-2631	4X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.10
	G380-2632	10X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.25



	G380-2633	40X INFINITY PLAN ACHROMAT OBJECTIVE, NA 0.65
	G380-2634	100X OIL IMMERSION INFINITY PLAN ACHROMAT OBJECTIVE, NA 1.25
	Miscellaneous Accessories	
	G380-5111	GRADUATED MECHANICAL STAGE (COAXIAL LOW POSITION CONTROLS) FOR G380 SERIES
	G380-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT FILTER HOLDER FOR G300 SERIES
	G380-7001	LED ILLUMINATOR PACK, SNAP-IN, INCLUDED 3 AA BATTERIES, FOR G380 SERIES ONLY
	G380-7002	MIRROR ASSEMBLY, SNAP-IN FOR G380 SERIES ONLY
	G380-7010	FIELD DIAGRAM KOEHLER SET FOR G380, G390 SERIES MICROSCOPES
	G380-8010	PHOTO LENS AND PHOTO ADAPTER WITH UNIVERSAL THREAD FOR G380 SERIES (REQ. CAMERA AND T-MOUNT)
	G380-8011	VERTICAL VIEWING ADAPTER (TO BE INSTALLED ON TRINOCULAR HEAD. REQUIRES EYEPIECE)
	G380-8031	PHASE KIT: INCLUDES 40XR SEMI-PLAN PHASE OBJECTIVE AND SLIDE-IN ANNULUS RING FOR G380 SERIES
	G380-8034	PHASE KIT: INCLUDES 10X SEMI-PLAN PHASE OBJECTIVE AND SLIDE-IN ANNULUS RING FOR G380 SERIES
	G380-8032	SLIDE-IN 40X PHASE ANNULUS RING FOR G380 SERIES
	G380-8907	HARD CARRYING CASE, PADDED WITH LOCK & KEY (2)
	G380-8001	DUST COVER FOR G380 SERIES
	G380-8906	DARKFIELD STOP FOR G380 SERIES
	G380-9008	BACKPACK™ FOR G380, G390 AND G490 SERIES MICROSCOPES, ACRYLIC ACCESORY ATTACHMENT
SPECIFICATION	Head	Trinocular
	Eyepiece	10X WF Hi-EP
	Objectives	Plan 4X/10X/40X/100X
	Condenser	NA1.25
	Diaphragm	Iris
	Stage	Mechanical



	Illumination	LED
	Bulb	3 Watt LED
	Focusing	Coaxial
SIZE AND WEIGHT	Shipping Size	17 in x 10.5 in x 17.5 in
	Weight	17 lbs (7.7 kg)
	Shipping Weight	18 lbs (8.2 kg)
WARRANTY		1 years

