

6 PLACE ROTOR (FOR FX/LX ONLY)

Model: C800-00



- ACCOMMODATES UP TO 6 SAMPLES
- BUILT FOR PRECISION AND RELIABLE PERFORMANCE
- DURABLE, HIGH-QUALITY CONSTRUCTION FOR LONG-TERM USE



A high-performance rotor designed for the FX and LX centrifuge models, offering space for 6 tubes. Ideal for small to medium sample batches.

ADDITIONAL INFORMATION

EQUIPMENT TYPE	Rotor
TOP SELECTION	No
SHIPPING WEIGHT	1 lb
ACCESSORIES	C800-02 Regular standard tube shield for 10ml -15ml tube, pack of 2 C800-03 Adapter for short tubes (75mm length, 4-5ml) C800-04 Adapter for pediatric or microtube (38-45mm tube length), pack of 2
WARRANTY	1 years

