

SHORT TUBE ADAPTERS/INSERTS, PK/2

Model: C800-6



- COMPATIBLE WITH 3-5ML BLOOD COLLECTION TUBES (75MM LENGTH)
- SET OF 2 ADAPTERS FOR CONVENIENCE
- PROVIDES A SECURE, LEAK-PROOF CONNECTION



This short tube adapter is designed to be used with 3-5ml blood collection tubes, specifically those with a 75mm length. Made of durable and high-quality materials, these adapters ensure a secure and stable connection, making them ideal for medical, laboratory, and clinical applications. The adapter fits seamlessly with blood collection tubes, providing an easy and efficient way to transfer or connect tubes for sample processing.



ADDITIONAL INFORMATION

EQUIPMENT TYPE	Adapter
MATERIAL	Wood
TOP SELECTION	No

