


DIGITAL POWERSPIN™ FX PRO

Model: C80508



- POWERSPIN™ FX PRO CENTRIFUGE
- TÜV RHEINLAND CERTIFIED
- FIXED SPEED 1000-3500 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES OR 16X15 ML TUBES



The C80508 PowerSpin™ FX PRO Centrifuge delivers consistent, fixed-speed operation at 3,400 RPM. Its 8-place rotor handles up to 8 x 10ml tubes, and the 99-minute digital timer ensures precise control. Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.



ADDITIONAL INFORMATION

TOP SELECTION	Yes	
ACCESSORIES	Accessories	
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
SPECIFICATION	Number of Tubes	8
	Tube Capacity	8-10 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Fixed
	Rotor	Fixed Angle
	Maximum RPM	3500
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	13 lbs
	Shipping Weight	15 lbs
WARRANTY	1 years	

