

# DIGITAL POWERSPIN™ PORTA-SPIN™ CENTRIFUGE

Model: C82608



- POWERSPIN™ PX CENTRIFUGE
- TÜV RHEINLAND CERTIFIED
- VARIABLE SPEED 1000-3500 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES



The PowerSpin™ Porta-Spin™ C82608 portable centrifuge is optimized for on-the-go applications. Running on 12VDC, it's perfect for fieldwork or locations lacking standard power outlets. The unit features a fixed 3,500 RPM, an 8-position rotor designed for 10mL tubes, and a 99-minute digital timer for precise control. Meeting high safety and electromagnetic compatibility standards, and certified by TÜV Rheinland, this centrifuge provides dependable performance in diverse environments.



**Every Purchase Includes a FREE  
AC Adapter – No Extra Cost!**



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	Yes
<b>ACCESSORIES</b>	<p><b>Accessories</b></p> <p>Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-75</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p>
<b>SPECIFICATION</b>	<p><b>Number of Tubes</b> 8</p> <p><b>Tube Capacity</b> 8-10 ml</p> <p><b>Blood</b> Yes</p> <p><b>Specimens Fecal</b> Not recommended</p> <p><b>Urine</b> Not recommended</p> <p><b>Microhematocrit</b> Not recommended</p> <p><b>Speed</b> Fixed</p> <p><b>Rotor</b> Fixed Angle</p> <p><b>Maximum RPM</b> 3500</p> <p><b>Maximum RCF(g)</b> 1397</p>
<b>SIZE AND WEIGHT</b>	<p><b>Size</b> 10.2 in x 11.8 in x 7.9 in</p> <p><b>Shipping Size</b> 13.4 in x 14.8 in x 10.6 in</p> <p><b>Weight</b> 13 lbs</p> <p><b>Shipping Weight</b> 15 lbs</p>
<b>WARRANTY</b>	1 years

