

# PRO SERIES 16 PLACE TUBE ROCKER

Model: L-TTR250



- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



The L-TTR250 PRO SERIES Tube Rocker/Mixer, part of a line also including the L-TTR150 and L-TTR450, offers high-throughput mixing for laboratory applications. This model accommodates up to 16 tubes and is designed to handle diverse tube sizes up to 16x120mm. It operates at a variable speed of 10-30 RPM ( $\pm 2$  RPM), allowing for gentle and thorough mixing. The unit is equipped with a non-slip, reversible, and autoclavable tube pad made of durable silicon rubber. The entire Pro Series line is TÜV Rheinland certified, ensuring compliance with safety regulations for reliable operation in any laboratory environment.



**ADDITIONAL INFORMATION**

<b>TOP SELECTION</b>	No
<b>SPECIFICATION</b>	<b>Capacity</b> 16 tubes <b>Speed</b> Variable <b>Voltage</b> 110V
<b>WARRANTY</b>	1 years

