

PRO SERIES 30 PLACE TUBE ROCKER

Model: L-TTR450



- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (± 2 RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



The L-TTR450 PRO SERIES Tube Rocker/Mixer, part of a line also including the L-TTR150 and L-TTR250, offers high-throughput mixing for laboratory applications. This model accommodates up to 30 tubes, while the Pro Series line also offers models for 32, 16, or 10 tubes to suit varying needs. Designed to handle diverse tube sizes up to 16x120mm, it operates at a variable speed of 10-30 RPM (± 2 RPM), allowing for gentle and thorough mixing. Features include a pause button and 110V operation for user convenience. The unit is equipped with a non-slip, reversible, and autoclavable tube pad made of durable silicon rubber. The entire Pro Series line is TÜV Rheinland certified, ensuring compliance with safety regulations for reliable operation in any laboratory environment.



ADDITIONAL INFORMATION

TOP SELECTION	No
SPECIFICATION	Capacity 30 tubes Speed Variable Voltage 110V
WARRANTY	1 years

