


BX ROTOR SAFETY LOCKING NUT

Model: C880-12



KEY FEATURES:

- ENHANCED SAFETY: SECURES THE ROTOR FIRMLY TO THE CENTRIFUGE SHAFT, MINIMIZING THE RISK OF ACCIDENTS.
- EASY INSTALLATION: SIMPLE AND STRAIGHTFORWARD TO INSTALL.



This essential safety component is a replacement part for BX series centrifuges. The safety locking nut ensures secure rotor attachment, preventing accidental disengagement during operation.



ADDITIONAL INFORMATION

MATERIAL	Steel
TOP SELECTION	No
WARRANTY	1 years

