

DIGITAL POWERSPIN™ PXV PRO

Model: C82708



- POWERSPIN™ PXV PRO CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES



The PXV Centrifuge offers a robust and dependable solution for laboratory sample separations. This benchtop centrifuge is engineered to deliver high-performance capabilities, reaching a maximum rotation speed of 3800 RPM and generating a relative centrifugal force of 1647×g. Designed with versatility in mind, it accommodates 8 x 15mL tubes, catering to a range of standard laboratory procedures. User-friendly digital controls allow for precise programming and monitoring of centrifugation runs, with adjustable speed and time settings, including a continuous operation option. The PXV operates with a quiet and efficient brushless DC motor, contributing to a more comfortable laboratory environment. This centrifuge meets top safety and electromagnetic compatibility standards for reliable operation in laboratory environments. It complies with UL 61010-1:2012, CAN/CSA-C22.2 NO.61010-1-12, and IEC 61010-2-020:2016 for safety. It also adheres to FCC Part 15 Subpart B: 2017 and IEC60950-1:2006 for electromagnetic compatibility, ensuring safe, interference-free use in labs globally.



ADDITIONAL INFORMATION

TOP SELECTION	Yes
ACCESSORIES	<p>Accessories</p> <p>Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-75</p> <p>Regular Tube Shield, White, pack of 2</p> <p>C8000-02</p> <p>5 ml Adapter for short tubes, pack of 8</p> <p>C8000-03-8</p>
SPECIFICATION	<p>Number of Tubes 8</p> <p>Tube Capacity 8-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Yes</p> <p>Urine Yes</p> <p>Microhematocrit With Adapters</p> <p>Speed Variable</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3800</p> <p>Maximum RCF(g) 1647</p>
SIZE AND WEIGHT	<p>Shipping Size 20 in x 16 in x 18 in</p> <p>Weight 21 lbs</p> <p>Shipping Weight 22 lbs</p>
WARRANTY	1 years

