


# PowerSpin™ FX Centrifuge

Model: C806



- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



The POWERSPIN™ FX Centrifuge is a practical and efficient solution for your centrifugation needs. This 110V centrifuge features a 6-place rotor, accommodating a variety of sample tubes. Its fixed speed of 3400 RPM and 30-minute digital timer simplify operation, while its durable construction ensures reliable performance in your laboratory. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins™ are pre-drilled/tapped to accept optional Tube Holdsters™



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>ACCESSORIES</b>	<p><b>Accessories (Sold Separately)</b></p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p>
<b>SPECIFICATION</b>	<p><b>Number of Tubes</b> 6</p> <p><b>Tube Capacity</b> 2-15 ml</p> <p><b>Blood</b> Yes</p> <p><b>Specimens Fecal</b> Not recommended</p> <p><b>Urine</b> Not recommended</p> <p><b>Microhematocrit</b> Not recommended</p> <p><b>Speed</b> Fixed</p> <p><b>Rotor</b> Fixed Angle</p> <p><b>Maximum RPM</b> 3400</p> <p><b>Maximum RCF(g)</b> 1320</p>
<b>SIZE AND WEIGHT</b>	<p><b>Size</b> 10.2 in x 11.8 in x 7.9 in</p> <p><b>Shipping Size</b> 13.4 in x 14.8 in x 10.6 in</p> <p><b>Weight</b> 13 lbs</p> <p><b>Shipping Weight</b> 15 lbs</p>



**WARRANTY**

1 years

